

Be Safe and Secure with the Bodypoint® Shower Chair Belt

THE BODYPOINT® SHOWER BELT

Air and water pass through the Aeromesh™ fabric providing a complete wash for the client and a short drying time for the belt.

Fastens with side-release buckle for the ease of caregivers.



Adjustable hole positions allow for easy length adjustment.



Rust-proof frame clamps allows you to secure the shower chair belt anywhere on the frame—even on an upholstered seat rail.

INTRODUCING THE BODYPOINT SHOWER CHAIR BELT

Bathing can be a dangerous business! The bathroom is commonly regarded as the most dangerous room in the home. Bath safety is a major concern for people with limited mobility; however, most shower chairs manufactured today do not come equipped with a safety belt to prevent shower chair users from slipping or falling.

The Bodypoint shower chair belt with our proprietary Aeromesh™ fabric provides safety without compromising

the comfort and health of the user. Water and air easily pass through the Aeromesh fabric, allowing the client to be bathed completely and the belt to dry quickly after use.

Its universal shape enables the belt to be used in a variety of ways depending on the needs of the client. The flexible, padded belt can be used across the chest or abdomen to prevent forward or lateral falls. Additionally, fitting a belt across the knees or legs can prevent clients from slipping down out of the shower chair.



Universal shape allows for use as a chest, abdominal or leg belt, as well as a calf strap.

PART NO.	DESCRIPTION
BB101M-22-1	Med, 22mm clamps
BB101M-25-1	Med, 25mm clamps
BB101L-22-1	Lg, 22mm clamps
BB101L-25-1	Lg, 25mm clamps
BB101M-1	Med, without clamps
BB101L-1	Lg, without clamps



2221 46TH AVE., LACHINE, QUEBEC H8T 3C9
 T. (514) 636-7575 • F. 877-636-8944
1-800-265-9830
www.advancedhealthcare.ca

