

# Comfortex® Xtenda Care™

## Affordable Excellence

### *Basic*

- Replaceable Stratos shear-reducing top cover is Healthcare's most durable and effective barrier against contamination.
- A layer of HR (high resiliency) top foam for enhanced comfort and pressure redistribution.
- A supportive foundation prevents heavier patients from bottoming out.

### *Firm Edging*

- High density, firm edging provides added support for easier patient transfers.
- Firm Edging is bonded to the Core and Top foam for greater strength and durability.

### *Secure Edging*

- Ultra Firm Edging rises above the horizon of the sleep surface to direct patients away from the edge of the bed and an accidental fall.
- Soft Central exit area facilitates supported transfers and sitting at the edge of the bed.
- A more conventional mattress appearance protects patient dignity.



"Great Value for a low price! These mattresses hold up better than most of the other mattresses we've had. The Residents seem very comfortable on the XtendaCare." Dir.of Nursing, AL

**Comfortex Innovating Surfaces**

1680 Wilkie Drive, Winona, MN 55987

1.800.445.4007

[www.prevent-injuries.com](http://www.prevent-injuries.com)

# Xtenda Care™

by Comfortex®

Provides Pressure Reduction and Comfort  
while giving you the options tailored to your needs.  
(Stages 1-2)

**Stratos Top Cover:** Healthcare's most durable guard against shear forces and moisture vapor pass-through.

**HR Top Foam:** High Resiliency material ensures superior cushioning properties, comfort and durability.

**Full 6.5" thickness:** Prevents bottoming-out.

**Hidden Handles:** Sewn to the underside of the mattress to guard against entanglement when transferring in and out of bed.



All Comfortex covers are water proof, vapor proof and anti-bacterial to maximize infection control.



EXCEL-MÉDICAL INC.

WWW.EXCEL-MEDICAL.CA

DISTRIBUÉ AU QUEBEC  
TEL.: 514-990-5196  
1-888-990-5196  
TELE.: 514-879-3496

Weight capacity: 250-300 lbs.  
Warranty: Full 3 Years