

# Fauteuil pour prélèvement Phlebotomy Chair

Conçu spécifiquement pour les prélèvements sanguins, ce fauteuil simple et confortable facilitera votre travail quotidien.

Pour faciliter le meilleur positionnement du bras, le fauteuil est muni d'appuis-bras légèrement concaves et ajustables en multi-position.

Le fauteuil est d'une hauteur siège-sol de 44cm (17.5 po) facilitant l'entrée et la sortie du fauteuil. Le dossier et le repose-pied sont synchronisés à l'aide d'un mécanisme permettant ainsi d'allonger le patient selon le besoin.

Especially design for blood drawing, this comfortable and practical chair will simplify your daily work load.

The chair has slightly concave, multidirectional positioning armrest to facilitate the positioning of the patient arms.

The chair has a seat to floor height of 44cm (17.5 in).

The back and the leg rest are synchronised by a mechanism to achieve a recline position when needed.



## SPÉCIFICATION TECHNIQUE

Dossier inclinable : 25-65° avec le repose-pied à l'horizontal  
Appui-bras : Ajustable en hauteur  
Ajustable en multi-position  
Recouvrement : Vinyle, protection antibactérienne facile à nettoyer  
Capacité de poids : 150 kg (330lbs)  
La structure : En acier et recouverte d'une peinture d'époxy.

## TECHNICAL SPECIFICATION

Folding backrest: 25-65° with horizontal footrest.  
Armrest: Foldable and height adjustable  
Multi-positioning adjustment  
Upholstery: Vinyl fabric, antibacterial and easy to clean  
Patient weights: Up to 150kg (330lbs)  
Structure: Steel with epoxy coating

## Distribué au Canada par:

Excel-Médical  
Succ. Place du Commerce, C.P.4  
Brossard, Qc Canada J4W 3L5  
Tél.: 514-990-5196 / 888-990-5196  
Courriel / email: info.1@excel-medical.ca  
www.excel-medical.ca

