


## 100 series



### MODEL 140 Standard Features

ST6030FT5SF-WW	Shower Taxi
External Dimensions	42"H x 22"W x 23"D
SSOF-GREY	Deluxe Padded Open Front Seat
C7	7 QT Pail
T5	Standard 5" Tente
DA	Drop Arm
SF-WW	Sliding Footrest with wheels
Weight Cap.	400 lbs

### MODEL 155 Standard Features

SC6013D	Shower Chair, Deluxe Elongated Commode Seat
External Dimensions	40"H x 22"W x 23"D
VACS	Vacuum Formed Seat 18"
C7	7 QT Pail
W4	Standard 4" Wagner
DA	Drop Arm
SF-WW	Sliding Footrest with wheels
Weight Cap.	300 lbs

### Optional Upgrades

MODEL #	Safety Belt	Privacy Skirt	Safety Bar	Soft Seat Open Front Grey	Soft Seat Closed Front White	Push Handles	White Pipe	Almond Pipe	Tilt Seat
Code #	B	PS	SB	SSOF-G	SSCF	PH	W	AP	TS
140	✓		✓		✓	✓	✓	✓	✓
155	✓			✓	✓	✓	✓	✓	✓

### Castor Upgrades

(available only on Model 155)

CUW4T3	Standard 4" to 3" Tente Casters (2locking/2non locking)
CUW4TT3	Standard 4" to 3" Tente Casters (4 locking)
CUW4TT5	Standard 4" to 5" Tente Casters (4 locking)

## 100 series



MODEL 120	Standard Features
SC6013D	Shower Chair Deluxe Elongated Commode Seat
External Dimensions	38.5"H x 22"W x 22.75"D
VACS	Vacuum Formed Seat (almond) 18"
C7	7 QT Pail
CA	Curved Arms
W3	Standard to 3" Wagner
Weight Cap.	300 lbs

MODEL 150	Standard Features
SC6013D	Shower Chair Deluxe Elongated Commode Seat
External Dimensions	38.5"H x 22"W x 22.75"D
VACS	Vacuum Formed Seat (almond) 21"
C7	7 QT Pail
DA	Drop Arm
W3	Standard to 3" Wagner
Weight Cap.	300 lbs

### Optional Upgrades

MODEL #	Safety Belt	Privacy Skirt	Safety Bar	Soft Seat Open Front Grey	Soft Seat Closed Front White	Push Handles	White Pipe	Almond Pipe	Tilt Seat
Code #	B	PS	SB	SSOF-G	SSCF	PH	W	AP	TS
120	✓			✓	✓	✓	✓	✓	✓
150	✓		✓	✓	✓	✓	✓	✓	✓

### Castor Upgrades (available on all models)

CUW3T3	Standard 3" to 3" Tente Casters (2locking/2non locking)
CUW3T3	Standard 3" to 3" Tente Casters (4locking)
CUW3T5-AL	Standard 3" to 5" Tente Casters (4locking)

\*All measurements are approximate. Custom modifications are available. Call for details.

## 100 series



MODEL 100	Standard Features
SC6013S	Shower Chair, Standard Commode Seat
External Dimensions	38.5"H x 22"W x 22.75"D
CA	Curved Arms
C7	Pail and Rail
W3	Wagner 3"
Weight Cap.	300 lbs



MODEL 105	Standard Features
SC6013S	Shower Chair, Standard Commode Seat
External Dimensions	38.5"H x 22"W x 22.75"D
CA	Curved Arms
C7	Pail and Rail
W3	Wagner 3"
F	Standard Footrest
Weight Cap.	300 lbs

### Optional Upgrades

MODEL #	Safety Belt	Privacy Skirt	Safety Bar	Soft Seat Open Front Grey	Soft Seat Closed Front White	Push Handles	White Pipe	Almond Pipe	Tilt Seat
Code #	<b>B</b>	<b>PS</b>	<b>SB</b>	<b>SSOF-G</b>	<b>SSCF</b>	<b>PH</b>	<b>W</b>	<b>AP</b>	<b>TS</b>
100	✓					✓	✓	✓	✓
105	✓					✓	✓	✓	✓

### Castor Upgrades (available on all models)

CUW3T3	Standard 3" to 3" Tente Casters (2locking/2non locking)
CUW3T3-AL	Standard 3" to 3" Tente Casters (4locking)
CUW3T5-AL	Standard 3" to 5" Tente Casters (4locking)

\*All measurements are approximate. Custom modifications are available. Call for details.

## 100 series

## Bariatric Chairs



### MODEL 160 Standard Features

SC6013W	Shower Chair, Deluxe Elongated Commode Seat
External Dimensions	38.5"H x 25"W x 25.5"D
VACS	Vacuum Formed Seat 21"
C7	7 QT Pail
W5	Standard 5" Wagner
Weight Cap.	350 lbs

### MODEL 165 Standard Features

SC6013W	Shower Chair, Deluxe Elongated Commode Seat
External Dimensions	38.5"H x 25"W x 25.5"D
VACS	Vacuum Formed Seat 21"
C7	7 QT Pail
W5	Standard 5" Wagner
CA	Curved Arms
SF-WW	Sliding Footrest with wheels
Weight Cap.	350 lbs

### MODEL 170 Standard Features

SC6014X	Shower Chair, Deluxe Elongated Commode Seat
External Dimensions	39.75"H x 28"W x 25.5"D
VACS	Vacuum Formed Seat 24"
C7	7 QT Pail
W5	Standard 5" Wagner
Optional	Belt
Weight Cap.	400 lbs

### Optional Upgrades

MODEL #	Safety Belt	Privacy Skirt	Safety Bar	Soft Seat Open Front Grey	Soft Seat Closed Front White	Push Handles	White Pipe	Almond Pipe	Tilt Seat
Code #	B	PS	SB	SSOF-G	SSCF	PH	W	AP	TS
160	✓	✓	✓	✓	✓	✓	✓	✓	✓
165	✓	✓	✓	✓	✓	✓	✓	✓	✓
170	✓	✓	✓	✓	✓	✓	✓	✓	✓

### Castor Upgrades

(available on all models)

CUW5T3	Standard 5" to 3" Tente Casters (2locking/2non locking)
CUW5TT3	Standard 5" to 3" Tente Casters (4locking)
CUW5TT5	Wagner 5" to 5" Tente Casters (4locking)

\*All measurements are approximate. Custom modifications are available. Call for details.