

Help Deter Resident Tampering



Deluxe Bed Sentry™ Model 91621

Features:

- Helps deter resident tampering with:
 - Dual button reset; must be pressed simultaneously to silence alarm
 - Volume and delay switches inside battery compartment
 - No on/off switch
- Patent pending wireless remote reset option
- Self-timing pad has visual indicator for pad expiration
- Built-in Sentry Pad cord strain relief
- Low battery indicator
- Uses single 9V battery
- Selectable 4-tone alarm
- Nurse call with remote alarm option
- DC adapter jack

Optional Accessories:

- Complete line of Sentry Pads
 - Bed or chair
 - Self-timing
 - Moisture-resistant
 - Damage-resistant
 - Seat belt
 - 45-day, 180-day, and 1-year options
- Integrated Latex-free cushion protector
- Universal bed bracket
- Universal chairstrap
- DC adapter
- Wireless remote reset pendant

Distributed by:

Excel-Médical Inc.

1.888.990.5196

Succ. Place du Commerce

Boite postale 4

Brossard, Québec J4W 3L5

www.excel-medical.ca

